

“My work requires me to hear conversations all day. Now I can hear speech more easily, even in large groups.”

“I am hearing things that I haven’t heard in years: the television, conversations, the voices of my grandchildren. It’s wonderful!”

“A lot of people did not even know that I was wearing a hearing aid. It is so small.”

A STYLE FOR EVERYONE

There’s a Conversa.NT style to suit every need and taste, from the tiniest completely-in-the-canal (CIC) to an expanded choice of behind-the-ear (BTE) models, including the ultra-small, stylish Moda. Your hearing healthcare professional will help you choose the best Conversa.NT for your needs.



ABOUT UNITRON HEARING

For 40 years, Unitron Hearing has been committed to making life better for people with hearing loss. To us, this means a commitment to not only develop high-quality hearing solutions, but to incorporate special features that really solve the everyday problems and concerns you have with hearing loss and hearing instruments. To you, this means a company that’s listening – and responding – to what matters in your life.

www.unitronhearing.com

conversa·NT™

GREAT CONVERSATIONS BEGIN HERE.





GREAT CONVERSATION PIECE

If you're looking for a hearing instrument that will optimize the way you hear conversations, even in noisy situations, Conversa.NT is right for you. Conversa.NT offers digital features that make everyday use simple and comfortable, and comes in a range of styles, sizes and colors to suit every taste.



DIFFERENTIATES SPEECH FROM NOISE

Listening to conversations in a restaurant can be challenging. Listening at the family dinner table can be just as challenging. Conversa.NT can help.

By constantly reading the sounds around you, Conversa.NT distinguishes the sounds of conversation from noise. When speech is present, speech enhancement engages to improve your ability to hear conversations. In noisy situations, noise reduction minimizes unwanted sounds and improves listening comfort.

HEAR WHAT YOU WANT

Conversa.NT keeps the focus on conversation, even when there's background noise. The directional microphones put priority on sounds directly in front of you (such as a person speaking to you), while detecting and suppressing noises from other directions, even if they're moving (such as a car driving by).

COMFORTABLE LISTENING

With Conversa.NT, your conversations won't be disrupted by annoying whistling. When something comes close to your ear (such as a hat or phone) you may experience feedback from your hearing aids. Conversa.NT's feedback canceller automatically controls feedback for more comfortable listening.

CONVERSA.NT FEATURES

Up to three listening programs – choose the best program for any environment

Directional microphones – emphasizes speech in front of you in noisy situations

Speech enhancement – boosts speech

Noise reduction – suppresses unwanted noise

Feedback canceller – reduces annoying whistling

Wind noise manager – provides comfort in windy environments

Ideal volume indicator – lets you know when you've reached your ideal volume