

unison™

The digital
choice



unitron
hearing 



It's all about
choice



You and Unison

Everyone loves choice, and when it comes to the Unison™ series of digital hearing instruments, the choice is all about you – your lifestyle, your listening preferences and your budget. The Unison series gives you the choice of three unique models, each offering the advantages of a high-quality digital hearing instrument, plus unique features to ensure the perfect choice for your needs.

Why Unison?

- ▶ High-value features and options work together to improve listening in more environments
- ▶ Clear, comfortable, more natural sound
- ▶ Choice of models to suit your hearing loss and lifestyle
- ▶ Wide choice of colors and sizes, including a new, virtually invisible style

Unique models for unique needs

Three unique Unison models address individual listening and lifestyle needs.

Unison 6

An ideal choice if you have an active lifestyle and interact in a wide variety of environments.

Unison 3

An ideal choice if you're exposed to fewer listening environments.

Unison Essential

An ideal choice if you're concerned about budget and/or you are making the transition from an analog to a digital hearing instrument.

The following pages will help you understand the differences between Unison 6, 3 and Essential. This comparison, along with the guidance of your hearing healthcare professional, will allow you to determine which Unison model is the best choice for you.



Choose your Unison™

Unison 6	Unison 3	Unison Essential
-------------	-------------	---------------------

Sound processing

What it is: during fitting, your hearing healthcare professional uses different sound processing strategies to address your specific hearing loss and your typical listening environments.

Benefit: the more sound processing strategies available, the more your hearing instruments can be customized to your hearing loss and listening preferences.

4 Sound strategies	2 Sound strategies	1 Sound strategy
--------------------------	--------------------------	------------------------

Feedback management

What it is: Unison's adaptive feedback canceller controls whistling, working continually to manage feedback in everyday situations as they occur.

Benefit: improved comfort and no concerns about embarrassing whistling during such everyday occurrences as pulling a hat over your ears or hugging.

•	•
---	---

Noise reduction

What it is: noise reduction analyzes your environment and diminishes background noise.

Benefit: improved listening comfort and decreased strain in noisy environments, as unwanted background sounds, like the hum of a dishwasher, are reduced.

•

	Unison 6	Unison 3	Unison Essential
<p>Directional microphones</p> <p>What it is: directional microphones put the emphasis on sounds directly in front of you rather than those beside or behind you.</p> <p>Benefit: because the sounds you want to hear, like conversations, are typically in front of you, directional microphones make it easier to focus on these sounds, even in noisy situations.</p>	•	•	
<p>Listening programs</p> <p>What it is: during fitting, your hearing healthcare professional can customize up to three listening situations based on your different listening environments and preferences.</p> <p>Benefit: as you move between listening environments or situations, simply push the program button to choose the best listening program for each situation. Unison confirms program changes with a discreet beep.</p>	3	3	1

Everyday conveniences

Everything about Unison is designed to make your life easier and your listening more enjoyable. These practical features and options are available in all three Unison models.



Ideal volume indicator – with a quiet beep, Unison lets you know precisely where your hearing healthcare professional has recommended you set your volume.

Low battery warning – when your batteries become low, Unison emits a discreet beep to give you advanced warning so you can replace them at your convenience.

Easy t-coil* – when you place the phone on your ear, easy t-coil automatically shifts to your telephone program. When you hang up, you'll be returned to the listening program previously in use.

Start-up mute – when you put your hearing instrument on, it won't begin amplifying sound until you push the program button. This prevents any whistling or feedback that might occur during start-up.

** Available as an option on all Unison in-the-ear styles. Telephone must be hearing instrument-compatible.*

Different styles for different needs

Unison comes in a wide range of discreet sizes and an array of colors. Based on your hearing loss and preferences, your hearing healthcare professional will help you select the right style for your needs.



Smaller size, big on comfort

The newest behind-the-ear (BTE) style of Unison is so small, lightweight and discreet it's virtually invisible. It fits comfortably behind your ear and uses thin, transparent and comfortable tubes that help keep your ear canal feeling more open. This style is best suited to certain types of hearing loss: your hearing healthcare professional can advise if it's right for you.







For 40 years, Unitron Hearing has been committed to making life better for people with hearing loss. To us, this means a commitment to not only develop high-quality hearing solutions, but to incorporate special features that really solve the everyday problems and concerns you have with hearing loss and hearing instruments. To you, this means a company that's listening – and responding – to what matters in your life.

Make the right choice with Unison.